

# Aerosol Containment Tent (ACT)

## Portable Negative Pressure

- Increase safety during aerosol generating procedures - intubation, extubation, BiPAP, CPAP, High Flow O<sub>2</sub>, and BVM Ventilation
- Increase negative pressure isolation capacity
- Safely expand non-invasive ventilation usage to avoid ICU admissions
- Enhanced flexibility in patient placement and temporary facilities
- Minimize risk of infectious aerosol transmission during elective procedures
- Provide continuity of patient care by avoiding procedure delays and cancellations



## Product Specifications


- 4 access ports for healthcare provider access
- Minimum -7.5 Pa pressure differential
- 99.97% efficiency HEPA filter
- Single-use disposal tent and duct
- Reusable frame, and fan filter unit
- Portable fan/filter on wheels
- Fits standard hospital beds and OR tables



MACH32


# Portable Negative Pressure

## Aerosols inside the ACT during simulated AGMPs



Aerosol Containment Tent (ACT)

## Aerosols at the physician during simulated AGMPs



> 7.5 Pa negative pressure

3000-4000 Air changes/hour

exceeds CSA guidelines for negative pressure isolation

## Ordering Information

Part Number	Description
MACH-4007, ACT V2 Starter Kit:	1 Fan/Filter Unit, 1 Reusable Frame, 25 Single-Use Tents
MACH-4001	25 Single-Use Tents
MACH-4004	25 Single-Use Ducts
Parts are available individually, please contact Andau Medical for additional information	



Andau Medical  
[sales@andaumedical.com](mailto:sales@andaumedical.com)  
 1-844-ANDAUM88  
[andaumedical.com](http://andaumedical.com)