

Medical Device Batteries (Health Canada Approved)

Device Model	Spec
Infusion Pump Batteries	
ALARIS MEDICAL SYSTEMS INC - 8015, 8110, 8220, 8300, 8600	12V 4AH NIMH PRESS CONTACTS
B BRAUN - 100ES, 200ES, 300ES, 400ES	12V 4.3AH NIMH .187 FASTON
MEDEX INC (Smiths Medical) - G-6000244, PROTEGE SYRINGE INFUSION PUMP - NIMH 1.7AH 3010, SYRINGE PUMP 3500, G-6001392-REV 2	7.2V 2.45AH NIMH PLUG
ABBOTT MEDICAL ELECTRONICS - 190706, LIFE CARE 175, 4100, 4P, 5500, 6000, PCA, XL Plum, XLD. BAXTER HEALTHCARE - Flo-gard 2000, 2100, 8000, 8500. HOSPIRA - H3B, PCA, PCA 3	8V 3.2AH SLA W/PLUG
B BRAUN - KANGAROO E PUMP	4.8V 3.6AH NIMH PLUG
B BRAUN - BASIC INFUSION PUMP VISTA	7.2V 2.2AH NICD W/ PLUG
BAXTER HEALTHCARE - 40A, AS50, AS50A, AS60A	6V 2100MAH NIMH
HOSPIRA - 830-03650-904	8V 3.2AH SLA PLUG
SMITHS MEDICAL (formerly MEDEX) & MEDFUSION SYSTEMS - SYRINGE INFUSOR 2001, 2010	7.2V 2000mAh NiCd
BAXTER HEALTHCARE - 40A, AS40, AS50, AS50A, AS60A	6V 1AH NICD PLUG
MEDRAD INC - SOLARIS EP MR INJECTION SYSTEM	12V 8AH SLA PURE .250 FASTON
Outlook 400ES,300ES,200ES,100ES, MED0169, R&D 6320, Unipower B11678	12V 4.5AH NIMH W/.187 FASTONS
Alphasource AS36296-I, MK BMED11676, R&D 6296-I, Spectrum 35083, spectrum 3570009, 35700bax2	7.4V 1.8AH 13WH LION INSERT
Monitor Batteries	
Welch Allyn 300 Series Vital Signs Monitor 53NTP	6V 4.5AH SLA W/PLUG
Welch Allyn LXI Spot (4500-84)	6V 6AH SLA W/ PLUG
Marquette Medical MAC1000, 1200	18V 2AH NICD PLUG
Spacelabs Medical - Monitor 1030, 1050, 90351, 90385	12V 650MAH NIMH PLUG
RESPIRONICS (aka HEALTHDYNE) - 418, 460009-027	6V 2.1AH NIMH PLUG
PROTOCOL SYSTEMS INC: PROPAQ MONITOR 200, 501-0008-00, 501-0008-00	
ENCORE 101EL, 102EL, 103EL, 104EL, 105EL, 106EL, 202EL, 204EL, 206EL, 244EL	8V 6.4AH SLA PLUG
Welch Allyn LXI Spot (4500-84) (Panasonic Cell)	6V 4.5AH SLA PLUG
CAS MEDICAL - NIBP 740-1, 740-2, 740-2T, 740-3	
VITAL SIGNS MONITOR 740, 750	7.2V 3.6AH NIMH PLUG
RESPIRONICS (aka HEALTHDYNE) - SMART MONITOR 900, 900S-10, 970SE-10	
OXIMETER 930, BATTERY PACK 915	6V 2.5AH NICD W/ PLUG
Welch Allyn LXI Spot (4500-84) (Panasonic Cell)	6V 4.5AH SLA PLUG PANASONIC
MEDIANA - YM1000, M6021-0	8.4V 7.2AH NIMH PLUG
Welch Allyn - 5200-84, BP/SPO2 5000, LIFESIGN ANALYZER 4200, 5000, 52000, SERIES LIFE SIGNS MONITOR 52000	6V 4.5AH SLA PLUG
8388, Dash Dash Series 2016989-002, -003, Flashpad, DR203, GE HealthCare Protransport, Stryker Prime NC2040XD, stretched integrated scale, Unipower B11823	10.8V 2.05AH 22WH LION
CardioPerfect - Series Electro Cardiograph CP150, Connex 65MTPX-B, Batt99, vital signs monitor 6000,6400,6500, R&D 6408,6408-a, Unipoer B11694	10.8V 6.75AH 73WH LION BATT99
Alphasource AS30019-1, Any battery SM201-6, DASH Series 20178757-002, 3000,4000,5000, DR Series DR201, GENEFORCE 2017857-002, MK BMED11194, R&D 5071,5071B, unipower B11194	11.1V 4.2AH 47WH LION DR201
Accutorr Series Monitors Plus,5 ; Alphasource As30278; Datascope Accutorr Series 0998-00-0444j81, Accutorr 5, Accutorr Plus Lion; Datascope Corp 0146-00-0069; Mindray 0146-00-0099, 115-018011-00, Pm- 8000, Pm-9000; Mk Battery Bmed11660; R&D Batteries 6132, 6132-A; Unipower B11641	11.1V 4.4AH 49WH LION MNDRY
VENTILATOR BATTERIES	
68339 AVEA VENT, B11407- UNIPOWER, VIASYS 68339K, VIASYS AVEA	24V 4.3AH NIMH PLUG
VIASYS 21542, VIASYS VELA	12V 4500MAH NIMH W/PLUG
Puritan-Nennet 840 vent, PB840 Ventilator,	24V 8.0AH SLA W/PLUG
Respironics (AKA HealthDYNE) 8-500016-00, BIPAP FOCUS VENTILATOR	9.6V 3.6AH NIMH PLUG
Nelcor Puritan-Bennett Achevia ventilator, Puritan Bennet L-007762-000	24V 2.3AH SLA .187 FASTON
Alphasource AS36212, Carefusion (19222-001, Ltv-1200,ltv-1000, ltv-900,ltv-950), Pulmonetice Systems (18608-001 Rev A, 10140-EP REVH), R& D 6212.	12V 4.5AH SLA PLUG
Coviden PB840, MED0024 (CR19058) , Puritan Bennet 840 VENT, PB840	24V 7.0AH SLA W/PLUG FR - GENESIS
Newport Medical instruments HT50 Ventilator	12V 2.2AH NIMH SOLDER TABS
Impact Medical 704-0754-01, 326 Portable Aspirator, 326M, UNIVENT EAGLE 754M,765,754,750	12V 4.5AH SLA PLUG
Siemens 300 Servo ventilator, 390 senoscreen, 9000C, 390 flatscreen, Is14500,servo300a,servo l memory,servo300memory,servo300 ventilator	3.6V 2.4AH LIT SIEMENS-MEMORY
Carefusion LTV-1200 Ventilator Internal battery, Pulmonetic Systems 18608-001 REV A, R& D 6212	12V 4.5AH SLA LTV 1200