

NARCOTICSCABINET


Management system for controlled substances and narcotics with a secure cabinet and integrated web-based inventory management software



✔ Highly secured cabinet with safe personal access



✔ Control of the entire drug issuance process



✔ Interface with hospital ERP and EHR systems

Benefits for hospitals

- Strict compliance with drug control regulations
- Individual drug pull-out and return
- Eliminates manual documentation
- Placed inside core areas and pharmacies
- Protects against theft and misuse
- Ability to track staff activities

Benefits for pharmacies

- Centralizes inventory management
- Flexible reporting system with alerts on expiring drugs and stockouts
- Precise, real-time consumption tracking
- Drug surplus returns and digital documentation destruction
- Documentation of broken ampoules or pull-out of many drugs at once



Our NARCOTICSCABINET solutions



ULTRA CABINET

Typically installed in OR core and can support up to 4 cabinets in one cluster. Up to 1,000 compartments in one cabinet and 67 "returned" compartments.



TABLETOP CABINET

For remote and specialized ORs, such as delivery room, ophthalmology ER, etc. Includes up to 200 compartments, 33 "returned" compartments.



Safe and simple personal access



Face Recognition



Biometric Identification




Worker Card



Personal Password


NARCOTICS CABINET | How it works



- 01** Physicians and other medical personnel are recognized with biometric, card or password identification.
- 02** Individual medications are dispensed for each patient by reading the patient barcode.
- 03** After dispensing, physicians report actual patient usage in the hospital clinical system.
- 04** Leftover items are returned to the cabinet; open dosages are eliminated and reported to the system.
- 05** The system alerts the pharmacist when inventory is running low by sending real-time notifications.
- 06** Pharmacists easily restock medications in the IDENTI secured cabinet.

Generate insights using IDENTIPLATFORM™ web server

- ▶ Centralized database for drug identification supported by advanced AI, machine-learning algorithms
- ▶ All cabinets are connected by a single network; selective access to reports according to privacy level
- ▶ Complete, real-time data visibility for precise management of all types of inventory



About IDENTI

IDENTI provides hospitals, medical device manufacturers and international logistics companies with powerful data-capture solutions to create a real-time account of medical inventory and consumption. Its multiple patent-protected technologies revolutionize hospital efficiency, increase patient safety, and lighten the burden on medical staff, ushering in a new user-friendly reality.

Other solutions

SNAP&GO™

Delivering visibility at point-of-use

TOTALSENSE

Implant and tissue tracking

KANBAN&PAR

Medical supplies replenishment

OPERATIONAL TOOLS & APPS

Track curtains and gas cylinders

