

# Airtraq

GUIDED VIDEO INTUBATION



Simple, effective,  
affordable video  
laryngoscopy

1.844.ANDAU88

[andaumedical.com](http://andaumedical.com)

[acutecare@andaumedical.com](mailto:acutecare@andaumedical.com)



ANDAU MEDICAL

# Airtraq sp

Adults  
Pediatrics  
Double Lumen  
Nasal (no channel)



Adult regular  
7.0-8.5 ETT  
M.O.\*: 16 mm  
**LDA011**



Adult small  
6.0-7.5 ETT  
M.O.\*: 15 mm  
**LDA021**



Double lumen  
28-41 Fr  
M.O.\*: 18 mm  
**LDA071**



Nasal  
M.O.\*: 15 mm  
**LDA061**

**NOW**  
4X the Image  
Resolution

Wi-Fi  
Camera



iOS / Android  
Mac / Windows



WiFi



Pediatric  
4.0-5.5 ETT  
M.O.\*: 12 mm  
**LDA031**



Infant  
2.5-3.5 ETT  
M.O.\*: 11 mm  
**LDA041**  
\* Mouth opening

**Fully Disposable  
One piece.**

**Airtraq improves routine intubations  
and it is specially effective for:**

Ready to use. VL for  
remote locations and  
emergencies.

- Cervical immobilization
- Obese
- Awake

- In combination with a fiberoptic
- Frontal or lateral intubations
- Poor dentition

# How Airtraq simplifies Video Intubation

- Airtraq's channel simplifies ETT insertion reducing intubation time.
- No stylet needed. Avoids potential injuries. Eliminates stylet costs.
- 90° shape minimizes hyperextension and reduces force required.
- Reduces hemodynamic reactivity.
- Facilitates intubation from any position (face to face, etc.).
- High intubation success rate in difficult airways.
- Easier to learn than other video laryngoscopes.

**Multiple  
viewing  
options**

**90° channel  
simplifies  
intubation**

**Minimum capital  
investment  
allows availability  
in every location**

# Visualization Options

for any Airtraq model



Wi-Fi  
Cam

Endo  
Cams



Direct  
View



Exclusively Distributed by:



[www.airtraq.com](http://www.airtraq.com)

Airtraq LLC  
2414 Lawton Ln.  
Rowlet, TX 75089, USA  
Phone: 1-877-624-7929  
E-mail: [info.usa@airtraq.com](mailto:info.usa@airtraq.com)



ANDAU MEDICAL

1.844.ANDAU88  
[andaumedical.com](http://andaumedical.com)  
[acutecare@andaumedical.com](mailto:acutecare@andaumedical.com)

Prodol Meditec S.A  
Muelle T. Olavarri, 5  
48930 Las Arenas, Spain  
Phone: +34 944 804 690  
E-mail: [info@airtraq.com](mailto:info@airtraq.com)